

LABORATORY STANDARD OPERATING PROCEDURE FOR PULSENET MLVA OF SHIGA TOXIN-PRODUCING <i>ESCHERICHIA COLI</i> O157 (STEC O157) AND <i>SALMONELLA ENTERICA</i> SEROTYPES TYPHIMURIUM AND ENTERITIDIS– APPLIED BIOSYSTEMS GENETIC ANALYZER 3130 PLATFORM	CODE: PNL23				
	Effective Date: <table border="1" style="display: inline-table; vertical-align: middle;"> <tr> <td style="width: 20px; text-align: center;">02</td> <td style="width: 20px; text-align: center;">26</td> <td style="width: 20px; text-align: center;">14</td> </tr> </table>			02	26
02	26	14			

1. PURPOSE: to describe the standardized laboratory protocol for molecular subtyping of Shiga toxin-producing *Escherichia coli* O157 (STEC O157) and *Salmonella enterica* serotypes Typhimurium and Enteritidis.

2. SCOPE: to provide the PulseNet participants with a single protocol for performing MLVA of STEC O157 and *Salmonella* serotypes Typhimurium and Enteritidis, thus ensuring inter-laboratory comparability of the generated results.

3. DEFINITIONS:

- 3.1 MLVA:** Multiple-locus variable-number tandem repeat analysis
- 3.2 VNTR:** Variable-number tandem repeat
- 3.3 DNA:** Deoxyribonucleic acid
- 3.4 DNase:** Deoxyribonuclenase
- 3.5 PCR:** Polymerase chain reaction
- 3.6 HPLC:** High purity liquid chromatography
- 3.7 dNTP:** Deoxyribonucleotide triphosphate
- 3.8 CDC:** Centers for Disease Control and Prevention
- 3.9 SOP:** Standard Operating Procedure

4. RESPONSIBILITIES/PROCEDURE

4.1. Biosafety warning: STEC O157 and *Salmonella* serotypes Typhimurium and Enteritidis with an infectious dose as low as 100 cells are human pathogens capable of causing serious disease. Always use a minimum of Biosafety level 2 practices and extreme caution when transferring and handling strains of these serotypes. Work in a biological safety cabinet when handling large amounts of cells. Disinfect or dispose of all plastic ware and glassware that come in contact with the cultures in a safe manner.

4.2. Reagents, supplies and equipment needed for DNA template preparation

- 4.2.1 Trypticase soy agar with 5 % sheep blood (TSA-SB) or comparable media
- 4.2.2 1 µl inoculation loops
- 4.2.3 0.5 ml microcentrifuge tubes
- 4.2.4 DNase-free, molecular biology -grade water
- 4.2.5 Vortex
- 4.2.6 Boiling water bath or thermocycler/thermal block accommodating 0.5 ml tubes
- 4.2.7 Tabletop centrifuge for high rpm (up to 13,000-14,000 rpm) spinning
- 4.2.8 Pipets (200 µl) for aliquoting 100 µl of DNase-free, molecular biology-grade water
- 4.2.9 Filtered Sterile Pipet tips

4.3. Reagents, supplies and equipment needed for PCR

- 4.3.1 DNA templates from isolates (keep at -20°C or -80°C freezer for long term)
- 4.3.2 PCR primers (see appendix PNL23-1)
 - 4.3.2.1 Fluorescent-labeled forward primers
 - 4.3.2.1.1 HPLC-purified
 - 4.3.2.2 Unlabeled reverse primers

VERSION:	REPLACED BY:	AUTHORIZED BY:	Page 1 of 26
-----------------	---------------------	-----------------------	---------------------

